INSTINCT[®]


Andros_INST0405127_INST0405128_INST0405129_INST0405126_INST0405125_TDS_V2

Technical Data Sheet | Modular Furniture

www.instinctproducts.co.uk

Andros 300mm Wall Column Two Door Unit

INST0405127, INST0405128, INST0405129, INST0405126, INST0405125



Roper Rhodes

Brassmill Lane Trading Estate, Bath, BA1 3JF. Technical Support – Tel: 01225 303 900 Email: sales@roperrhodes.co.uk

All information on the data sheet is correct at the time of publication. All measurements are approximate and subject to standard tolerances.

INSTINCT[®]

Issue No. v2 07 2024

Andros_INST0405127_INST0405128_INST0405129_INST0405126_INST0405125_TDS_V2

Technical Data Sheet | Modular Furniture

www.instinctproducts.co.uk

Andros 300mm Wall Column Two Door Unit

INST0405127, INST0405128, INST0405129, INST0405126, INST0405125

FAQs

Q - Are the hinges soft close?
A - The hinges are both soft close and quick release.
Q - Is the unit reversable?

A – Yes, there are hangers on the top and bottom to reverse the unit. See instructions.

- Q Are the shelves removable?
- A No. There are two fixed shelves
- Q Can the doors be adjusted?
- A The doors of the unit can be adjusted as detailed in the supplied instructions.
- Q How deep is the unit when the doors are fully open?
- A Approx. 550mm including handle. (approx 90°).

Q - How do I fix the unit to the wall? A - There are adjustable wall hangers in the unit and wall plate kit is provided as detailed in the supplied instructions.

Roper Rhodes

Brassmill Lane Trading Estate, Bath, BA1 3JF. Technical Support – Tel: 01225 303 900 Email: sales@roperrhodes.co.uk

All information on the data sheet is correct at the time of publication. All measurements are approximate and subject to standard tolerances.